


## Generate spindle poles tracks from ImageJ FIJI (1.52v) plugin TrackMate


1. Film L4 worms as described (Zellag, 2021) with a confocal SD microscope.
2. Track the spindle poles as described (Zellag, 2021) or with the ImageJ plugin TrackMate as followed:
  - 2.1. Open your registered tif file in FIJI imageJ
  - 2.2. Launch TrackMate from the menu *Plugins* › *Tracking* › *TrackMate*




- 2.3. The settings are usually correct, click “Next”
- 2.4. Select log detector and click “Next”




- 2.5. Select the channel corresponding to tubulin, set the estimated blob diameter to 2,5  $\mu\text{m}$  and the threshold to 10. Make sure the “do subpixel localization” is checked.



- 2.6. Click Next
- 2.7. Let the detection process run.
- 2.8. Click Next




2.9. Set the initial filter to “auto” (note here that you could do it manually depending of the quality of the movie, in a manner to catch most of the spindle poles.



Click Next.


- 2.10. Select “Hyperstack displayed” and click Next
- 2.11. Do not set additional filter and click Next




2.12. Select “Simple LAP tracker” and

2.13. click Next


2.14. Set “Linking max distance” to 5.0 micros, “Gap-closing max distance” to 5.0



micron and “Gap-closing max frame gap” to 2 and click Next




2.15. Click Next again




2.16. Do not set additional filters on tracks and click Next




help of the movie. Each time you click on a point of this TrackScheme window, the point is highlighted on the movie and vice-versa.



2.19. Delete all the irrelevant tracks and make sure there is no empty spot in the track (make sure it is complete)




2.20. Click Next




2.21. Click Next






2.22. Select « Trim non-visible data » and click Execute



2.23. Click Save as movie\_name.xml eg: 2020-07-16\_UM776\_L4\_26C\_24h-01.xml



2.24. Select « Export all spots statistics » and

2.25. click Execute

2.26. Save all spots statistics as “movie\_name.txt (eg: 2020-07-16\_UM776\_L4\_26C\_24h-01.txt)

2.27. You can now close ImageJ.